

Anti-GPSM2 antibody (Internal) (STJ71084)

STJ71084

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-GPSM2 (Internal) is suitable for use in ELISA, Immunohistochemistry and Immunofluorescence research applications.
Applications	Pep-ELISA/IHC/IF
Host/Source	Goat
Reactivity	Human/Dog

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000. IHC: Paraffin embedded Human Testis. Recommended concentration: 5µg/ml. IF: Strong expression of the protein seen in the cytoplasm/centrosomes of A431 cells. Recommended concentration: 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Formulation	
Isotype	IgG
Storage	Store at -20°C on receipt and minimise freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	29899
Gene Symbol	GPSM2
Uniprot ID	GPSM2_HUMAN
Immunogen	
Immunogen	Internal
Region	
Specificity	
Immunogen	NRLKGKKYKTNSSTK
Sequence	