

Anti-CLPB antibody (C-Term) (STJ71168)

STJ71168

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-CLPB (C-Term) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave bands at approx 50kDa and 48kDa in lysates of HeLa and HepG2 after 0.1µg/ml antibody staining. Please note that currently we cannot find an explanation in
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81570
Gene Symbol	CLPB
Uniprot ID	CLPB_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	DKDSKTRRLDIRAP