

Anti-Transmembrane protein 37 antibody (Internal) (STJ71308)

STJ71308

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Transmembrane protein 37 (Internal) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave an approx 26-28kDa band in Human Brain (Cerebellum, Frontal Cortex and Hippocampus) lysates after 0.1µg/ml antibody staining. Please note that currently we
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	140738
Gene Symbol	TMEM37
Uniprot ID	CCGL_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	TNQTCFRDLGQ