

Anti-TRPV5 antibody (Internal) (STJ71557)

STJ71557

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-TRPV5 (Internal) is suitable for use in ELISA, Western Blot, Immunocytochemistry and Immunofluorescence research applications.
Applications	Pep-ELISA/WB/ICC/IF
Host/Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Pig/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:16000. WB: Approx 85kDa band observed in Human Pancreas lysates (calculated MW of 82.6kDa according to NP_062815.2). Recommended concentration: 1-3µg/ml. Primary incubation was 1 hour. This antib
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	56302
Gene Symbol	TRPV5
Uniprot ID	TRPV5_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	SHRGWEILRQNT



STJ71557 (1µg/ml) staining of Human Pancreas lysate (35µg protein in RIPA buffer). Detected by chemiluminescence.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081