

Anti-APOLD1/VERGE antibody (C-Term) (STJ71655)

STJ71655

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-APOLD1/VERGE (C-Term) is suitable for use in ELISA, Western Blot and Flow Cytometry research applications.
Applications	Pep-ELISA/WB/FC
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000. WB: Approx. 35kDa band observed in lysates of cell lines A431, Jurkat and U251 (calculated MW of 30.5kDa according to Human NP_001123887.1). This molecular weight is routinely observed by ot
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81575
Gene Symbol	APOLD1
Uniprot ID	APLD1_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	KSSRGHDLKISADQ