

Anti-EBI3 antibody (Internal) (STJ71763)

STJ71763

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-EBI3 (Internal) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave an approx 85kDa band in lysates of cell lines Daudi, Jurkat and MOLT4 after 0.02µg/ml antibody staining. Please note that currently we cannot find an explanation.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10148
Gene Symbol	EBI3
Uniprot ID	IL27B_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	IRYKRQGAARFHR