

Anti-SAFB2 antibody (Internal) (STJ71869)

STJ71869

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-SAFB2 (Internal) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Preliminary experiments gave bands at approx 75kDa and 60kDa in lysates of HeLa cell lines after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in t
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9667
Gene Symbol	SAFB2
Uniprot ID	SAFB2_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	This antibody is not expected to cross-react with SAFB1.
Immunogen Sequence	GEHNQRLEEHQARA