

Anti-DTX1/Deltex-1 antibody (Internal) (STJ72066)

STJ72066

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-DTX1/Deltex-1 (Internal) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: Preliminary experiments gave an approx 48kDa band in Mouse Fetal Brain lysates and in lysates of NIH3T3 after 1µg/ml antibody staining. Please note that currently we cannot find an expl
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1840
Gene Symbol	DTX1
Uniprot ID	DTX1_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	CHHIENVLKEDAR