

Anti-ZDHHC1 antibody (174-182) (STJ72387)

STJ72387

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-ZDHHC1 (174-182) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:2000. WB: Preliminary experiments gave bands at approx 48kDa and 20kDa in lysates of cell line HEK293 after 2µg/ml antibody staining. Please note that currently we cannot find an explanation in t
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29800
Gene Symbol	ZDHHC1
Uniprot ID	ZDHC1_HUMAN
Immunogen	
Immunogen Region	174-182
Specificity	
Immunogen Sequence	RNYRLFLHS