

Anti-BTBD7 antibody (N-Term) (STJ72469)

STJ72469

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-BTBD7 (N-Term) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA/WB
Host/Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Pig/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:64000.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55727
Gene Symbol	BTBD7
Uniprot ID	BTBD7_HUMAN
Immunogen	
Immunogen Region	N-Term
Specificity	This antibody is expected to recognize both reported isoforms (NP_001002860.2; NP_060637.1).
Immunogen Sequence	ESKLYSLDHGHEKPQ



MDCK overexpressing Mouse Btbd7 (third lane) and probed with STJ72469 (mock transfection in first lane, full construct chemically suppressed in second lane).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081