

Anti-Trib1-Mouse antibody (304-317) (STJ72655)

STJ72655

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-Trib1-Mouse (304-317) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA/WB
Host/Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:8000. WB: Approx 42kDa band observed in lysates of cell line NIH3T3 (calculated MW of 41.3kDa according to NP_653132.1). Recommended concentration: 0.1-0.3µg/ml.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10221
Gene Symbol	TRIB1
Uniprot ID	TRIB1_HUMAN
Immunogen	
Immunogen Region	304-317
Specificity	
Immunogen Sequence	PEHVSPKARCLIRS



STJ72655 (0.1 µg/ml) staining of NIH3T3 lysate (35 µg protein in RPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.