

Anti-POLRMT antibody (C-Term) (STJ73010)

STJ73010

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-POLRMT (C-Term) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Preliminary experiments gave bands at approx 90kDa and 35kDa in Human HepG2 lysates after 1µg/ml antibody staining. Please note that currently we cannot find an explanation in the lit
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5442
Gene Symbol	POLRMT
Uniprot ID	RPOM_HUMAN
Immunogen	
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	KETLQAVPKPGA