

Anti-CD3EAP antibody (Internal) (STJ73099)

STJ73099

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-CD3EAP (Internal) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:4000. WB: Preliminary experiments in nuclear lysates of cell lines HeLa, Jurkat and A431 gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate a
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10849
Gene Symbol	POLR1G
Uniprot ID	RPA34_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	PEDKTVKQEINTEP