

Anti-LRP1 antibody (3941-3951) (STJ73142)

STJ73142

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-LRP1 (3941-3951) is suitable for use in ELISA research applications.
Applications	Pep-ELISA
Host/Source	Goat
Reactivity	Human/Mouse/Rat/Dog/Pig/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:128000. WB: Not yet tested-our routinely used western blotting protocol does not allow detection of proteins as large as the calculated size of 505kDa according to NP_002323.2. Therefore we cannot
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4035
Gene Symbol	LRP1
Uniprot ID	LRP1_HUMAN
Immunogen	
Immunogen Region	3941-3951
Specificity	
Immunogen Sequence	HRRQIDRGVTH