

Anti-HOXA10, Biotinylated antibody (Internal) {Biotin} (STJ73545)

STJ73545

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-HOXA10, Biotinylated (Internal) is suitable for use in ELISA and Western Blot research applications.
Applications	Pep-ELISA/WB
Host/Source	Goat
Reactivity	Human/Mouse/Pig/Cow

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Biotin
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Dilution Range	Peptide ELISA: antibody detection limit dilution 1:32000. WB: Approx 38kDa band observed in Human Skeletal Muscle lysates (calculated MW of 42.4kDa according to NP_061824.3). Recommended concentration: 0.2-0.6µg/ml.
Formulation	0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin. NA
Isotype	IgG
Storage Instruction	Store at -20°C on receipt and minimise freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3206
Gene Symbol	HOXA10
Uniprot ID	HXA10_HUMAN
Immunogen	
Immunogen Region	Internal
Specificity	
Immunogen Sequence	PQATSCSFAQNIKE.



Biotinylated STJ73033 (0.2µg/ml) staining of Human Skeletal Muscle lysate (35µg protein in RIPA buffer). Primary incubation was 1 hour. Detected by chemiluminescence.