

Anti-Phospho-CaMKII Alpha/Beta/Delta-Thr305 antibody (271-320 aa) (STJ90439)

STJ90439

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Calcium/calmodulin-dependent protein kinase type II subunit alpha and Calcium/calmodulin-dependent protein kinase type II subunit beta and Calcium/calmodulin-dependent protein kinase type II subunit delta-Thr30
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	IHC 1:100-1:300 ELISA 1:5000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	816 817 815 CAMK2B CAMK2D KCC2B_HUMAN KCC2D_HUMAN KCC2A_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human CaMK2 alpha/beta/delta around the phosphorylation site of Thr305 at the amino acid range 271-320
Immunogen Region	271-320 aa
Specificity	Phospho-CaMKII Alpha/Beta/Delta (T305) Polyclonal Antibody detects endogenous levels of CaMKII Alpha/Beta/Delta protein only when phosphorylated at T305.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081