

Anti-Phospho-E2F1-Thr433 antibody (388-437 aa) (STJ90458)

STJ90458

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Transcription factor E2F1-Thr433 (388-437 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1869
Gene Symbol	E2F1
Uniprot ID	E2F1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human E2F1 around the phosphorylation site of Thr433 at the amino acid range 388-437
Immunogen Region	388-437 aa
Specificity	Phospho-E2F-1 (T433) Polyclonal Antibody detects endogenous levels of E2F-1 protein only when phosphorylated at T433.
Immunogen Sequence	



Western blot analysis of lysates from HeLa cells treated with Etoposide 25uM/24h, using E2F1 (Phospho-Thr433) Antibody. The lane on the right is blocked with the phospho peptide.