

Anti-Phospho-GAB2-Tyr452 antibody (STJ90630)

STJ90630

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-GRB2-associated-binding protein 2-Tyr452 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	9846
Gene Symbol	GAB2
Uniprot ID	GAB2_HUMAN
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human Gab 2 (phospho Tyr452)
Immunogen Region	
Specificity	Phospho-Gab 2 (Y452) Polyclonal Antibody detects endogenous levels of Gab 2 protein only when phosphorylated at Y452.
Immunogen Sequence	