

Anti-ACAD10 antibody (231-280 aa) (STJ91433)

STJ91433

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acyl-CoA dehydrogenase family member 10 (231-280 aa) is suitable for use in Immunohistochemistry, Immunofluorescence, Western Blot and ELISA research applications.
Applications	IHC/IF/WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Range	IHC 1:100-1:300 ELISA 1:40000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	80724
Gene Symbol	ACAD10
Uniprot ID	ACD10_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human ACAD10 at the amino acid range 231-280
Immunogen Region	231-280 aa
Specificity	ACAD-10 Polyclonal Antibody detects endogenous levels of ACAD-10 protein.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human heart, using ACAD10 Antibody. The picture on the right is blocked with the synthesized peptide.