

Anti-AMPD2 antibody (131-180 aa) (STJ91577)

STJ91577

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AMP deaminase 2 (131-180 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	271
Gene Symbol	AMPD2
Uniprot ID	AMPD2_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human AMPD2 at the amino acid range 131-180
Immunogen Region	131-180 aa
Specificity	AMPD2 Polyclonal Antibody detects endogenous levels of AMPD2 protein.
Immunogen Sequence	



Western blot analysis of the lysates from HeLa cells using AMPD2 antibody.



Western blot analysis of various cells using AMPD2 Polyclonal Antibody

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081