

Anti-ARFGAP3 antibody (280-360 aa) (STJ91672)

STJ91672

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ADP-ribosylation factor GTPase-activating protein 3 (280-360 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	26286
Gene Symbol	ARFGAP3
Uniprot ID	ARFG3_HUMAN
Immunogen	Synthesized peptide derived from ARF GAP3 at the amino acid range 280-360
Immunogen Region	280-360 aa
Specificity	ARF GAP3 Polyclonal Antibody detects endogenous levels of ARF GAP3 protein.
Immunogen Sequence	