

Anti-ARHGAP12 antibody (310-390 aa) (STJ91680)

STJ91680

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Rho GTPase-activating protein 12 (310-390 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	94134
Gene Symbol	ARHGAP12
Uniprot ID	RHG12_HUMAN
Immunogen	Synthesized peptide derived from ARHGAP12 at the amino acid range 310-390
Immunogen Region	310-390 aa
Specificity	ARHGAP12 Polyclonal Antibody detects endogenous levels of ARHGAP12 protein.
Immunogen Sequence	