

Anti-Casein Kinase I Gamma 1/2/3 antibody (229-278 aa) (STJ92013)

STJ92013

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Casein kinase I isoform gamma-1 and Casein kinase I isoform gamma-2 and Casein kinase I isoform gamma-3 (229-278 aa) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC 1:100-1:300
Range	ELISA 1:40000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID [1456](#)
[1455](#)
[53944](#)
[CSNK1G3](#)
[CSNK1G2](#)
[KC1G3_HUMAN](#)
[KC1G2_HUMAN](#)
[KC1G1_HUMAN](#)

Immunogen The antiserum was produced against synthesized peptide derived from the human CK-1gamma1/2/3 at the amino acid range 229-278
Immunogen Region 229-278 aa
Specificity Casein Kinase I Gamma 1/2/3 Polyclonal Antibody detects endogenous levels of Casein Kinase I Gamma 1/2/3 protein.
Immunogen Sequence



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4A-C overnight). High-pressure and temperature Tris-EDTA, pH8.0 was used for antigen retrieval. Negative control (right) obtained from antibody was pre-absorbed by immunogen peptide.