

Anti-CD179b antibody (26-75 aa) (STJ92099)

STJ92099

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Immunoglobulin lambda-like polypeptide 1 and Immunoglobulin lambda constant 1 and OBSOLETE and OBSOLETE and Immunoglobulin lambda constant 6 and Immunoglobulin lambda constant 7 (26-75 aa) is suitable for use in
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3543 NA NA NA 3543 NA NA NA
Gene Symbol	IGLL1 IGLC6 IGLL1_HUMAN IGLC6_HUMAN IGLC1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human CD179b at the amino acid range 26-75
Immunogen Region	26-75 aa
Specificity	CD179b Polyclonal Antibody detects endogenous levels of CD179b protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081