

Anti-CELSR3 antibody (3195-3244 aa) (STJ92216)

STJ92216

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cadherin EGF LAG seven-pass G-type receptor 3 (3195-3244 aa) is suitable for use in Immunofluorescence and ELISA research applications.
Applications	IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IF 1:200-1:1000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1951
Gene Symbol	CELSR3
Uniprot ID	CEL3_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human CELSR3 at the amino acid range 3195-3244
Immunogen Region	3195-3244 aa
Specificity	CELSR3 Polyclonal Antibody detects endogenous levels of CELSR3 protein.
Immunogen Sequence	



Immunofluorescence analysis of A549 cells, using CELSR3 Antibody. The picture on the right is blocked with the synthesized peptide.