

Anti-CSGALNACT2 antibody (61-110 aa) (STJ92494)

STJ92494

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Chondroitin sulfate N-acetylgalactosaminyltransferase 2 (61-110 aa) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC 1:100-1:300
Range	ELISA 1:40000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	55454
Gene Symbol	CSGALNACT2
Uniprot ID	CGAT2_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human CSGALNACT2 at the amino acid range 61-110
Immunogen Region	61-110 aa
Specificity	CsGalNACT-2 Polyclonal Antibody detects endogenous levels of CsGalNACT-2 protein.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using CSGALNACT2 Antibody. The picture on the right is blocked with the synthesized peptide.