

Anti-CYP26C1 antibody (Internal) (STJ92565)

STJ92565

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytochrome P450 26C1 (Internal) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	340665
Gene Symbol	CYP26C1
Uniprot ID	CP26C_HUMAN
Immunogen	Synthesized peptide derived from the Internal region of human CYP26C1.
Immunogen Region	Internal
Specificity	CYP26C1 Polyclonal Antibody detects endogenous levels of CYP26C1 protein.
Immunogen Sequence	