

Anti-DMGDH antibody (817-866 aa) (STJ92727)

STJ92727

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Dimethylglycine dehydrogenase, mitochondrial (817-866 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29958
Gene Symbol	DMGDH
Uniprot ID	M2GD_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human DMGDH at the amino acid range 817-866
Immunogen Region	817-866 aa
Specificity	DMGDH Polyclonal Antibody detects endogenous levels of DMGDH protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081