

Anti-EphB1/2/3/4 antibody (566-615 aa) (STJ92950)

STJ92950

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ephrin type-B receptor 1 and Ephrin type-B receptor 2 and Ephrin type-B receptor 3 and Ephrin type-B receptor 4 (566-615 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA resea
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID [2047](#)
[2048](#)
[2050](#)
[EPHB1](#)
[EPHB2](#)
[EPHB1_HUMAN](#)
[EPHB2_HUMAN](#)
[EPHB4_HUMAN](#)

Immunogen The antiserum was produced against synthesized peptide derived from the human EPHB1/2/3/4 at the amino acid range 566-615
Immunogen Region 566-615 aa
Specificity EphB1/2/3/4 Polyclonal Antibody detects endogenous levels of EphB1/2/3/4 protein.
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081