

## Anti-KHSRP antibody (N-Term) (STJ93050)

STJ93050

### GENERAL INFORMATION


|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Far upstream element-binding protein 2 (N-Term) is suitable for use in Western Blot and ELISA research applications. |
| <b>Applications</b>      | WB/ELISA   |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human/Mouse/Rat  |

### PRODUCT PROPERTIES

|                       |   |
|-----------------------|---|
| <b>Clonality</b>      | Polyclonal  |
| <b>Clone ID</b>       |   |
| <b>Concentration</b>  | 1 mg/mL   |
| <b>Conjugation</b>    | Unconjugated  |
| <b>Purification</b>   | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b> | WB 1:500-1:2000<br>ELISA 1:10000  |
| <b>Formulation</b>    | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.   |
| <b>Isotype</b>        | IgG   |
| <b>Storage</b>        | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |
| <b>Instruction</b>    |   |

### TARGET INFORMATION

|                           |   |
|---------------------------|---|
| <b>Gene ID</b>            | 8570  |
| <b>Gene Symbol</b>        | KHSRP   |
| <b>Uniprot ID</b>         | FUBP2_HUMAN   |
| <b>Immunogen</b>          | Synthesized peptide derived from the N-terminal region of human FBP2. |
| <b>Immunogen Region</b>   | N-Term  |
| <b>Specificity</b>        | FBP2 Polyclonal Antibody detects endogenous levels of FBP2 protein.   |
| <b>Immunogen Sequence</b> |   |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081