

Anti-FoxD4/D4L antibody (281-330 aa) (STJ93110)

STJ93110

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Forkhead box protein D4 and Forkhead box protein D4-like 4 and Forkhead box protein D4-like 5 and Forkhead box protein D4-like 6 and Forkhead box protein D4-like 3 and Forkhead box protein D4-like 6 (281-330 aa) is sui
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	286380 2298 349334
Gene Symbol	FOXD4L3 FOXD4 FX4L3_HUMAN FOXD4_HUMAN FX4L4_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human FOXD4/L2/L3/L4/L5/L6 at the amino acid range 281-330
Immunogen Region	281-330 aa
Specificity	FoxD4/D4L Polyclonal Antibody detects endogenous levels of FoxD4/D4L protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081