

Anti-GPT2 antibody (190-270 aa) (STJ93408)

STJ93408

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Alanine aminotransferase 2 (190-270 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84706
Gene Symbol	GPT2
Uniprot ID	ALAT2_HUMAN
Immunogen	Synthesized peptide derived from GPT2 at the amino acid range 190-270
Immunogen Region	190-270 aa
Specificity	GPT2 Polyclonal Antibody detects endogenous levels of GPT2 protein.
Immunogen Sequence	