

Anti-Histone H1 antibody (1-50 aa) (STJ93515)

STJ93515

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Histone H1.5 and Histone H1.3 and Histone H1.4 (1-50 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID [3007](#)
[3008](#)
[3009](#)
[H1-3](#)
[H1-4](#)
[H13_HUMAN](#)
[H14_HUMAN](#)
[H15_HUMAN](#)
[<](#)

Immunogen The antiserum was produced against synthesized peptide derived from the human Histone H1 at the amino acid range 1-50
Immunogen Region 1-50 aa
Specificity Histone H1 Polyclonal Antibody detects endogenous levels of Histone H1 protein.
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081