

Anti-IL11RA antibody (300-380 aa) (STJ93675)

STJ93675

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Interleukin-11 receptor subunit alpha (300-380 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3590
Gene Symbol	IL11RA
Uniprot ID	I11RA_HUMAN
Immunogen	Synthesized peptide derived from IL-11R Alpha at the amino acid range 300-380
Immunogen Region	300-380 aa
Specificity	IL-11R Alpha Polyclonal Antibody detects endogenous levels of IL-11R Alpha protein.
Immunogen Sequence	



Western blot analysis of various cells using IL-11R Alpha Polyclonal Antibody

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081