

Anti-LYL1 antibody (151-200 aa) (STJ93972)

STJ93972

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein lyl-1 (151-200 aa) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4066
Gene Symbol	LYL1
Uniprot ID	LYL1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human Lyl-1 at the amino acid range 151-200
Immunogen Region	151-200 aa
Specificity	Lyl-1 Polyclonal Antibody detects endogenous levels of Lyl-1 protein.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using Lyl-1 Antibody. The picture on the right is blocked with the synthesized peptide.