

Anti-NKX2-4 antibody (10-59) (STJ94503)

STJ94503

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Homeobox protein Nkx-2.4 (10-59) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	644524
Gene Symbol	NKX2-4
Uniprot ID	NKX24_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human NKX24. AA range:10-59
Immunogen Region	10-59
Specificity	Nkx-2.4 Polyclonal Antibody detects endogenous levels of Nkx-2.4 protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081