

Anti-NUCKS1 antibody (28-77) (STJ94572)

STJ94572

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Nuclear ubiquitous casein and cyclin-dependent kinase substrate 1 (28-77) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC 1:100-1:300
Range	ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	64710
Gene Symbol	NUCKS1
Uniprot ID	NUCKS_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human NUCKS1. AA range:28-77
Immunogen Region	28-77
Specificity	NUCKS Polyclonal Antibody detects endogenous levels of NUCKS protein.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma tissue using NUCKS1 Antibody. The picture on the right is blocked with the synthesized peptide.