

Anti-OAT antibody (100-180) (STJ94584)

STJ94584

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ornithine aminotransferase, mitochondrial (100-180) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4942
Gene Symbol	OAT
Uniprot ID	OAT_HUMAN
Immunogen	Synthesized peptide derived from OAT. at AA range: 100-180
Immunogen Region	100-180
Specificity	OAT Polyclonal Antibody detects endogenous levels of OAT protein.
Immunogen Sequence	