

Anti-POFUT2 antibody (361-410) (STJ94603)

STJ94603

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-GDP-fucose protein O-fucosyltransferase 2 (361-410) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23275
Gene Symbol	POFUT2
Uniprot ID	OFUT2_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human POFUT2. AA range:361-410
Immunogen Region	361-410
Specificity	O-FucT-2 Polyclonal Antibody detects endogenous levels of O-FucT-2 protein.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human placenta tissue, using POFUT2 Antibody. The picture on the right is blocked with the synthesized peptide.