

Anti-OR4F17 antibody (51-100) (STJ94709)

STJ94709

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Olfactory receptor 4F4/4F5/4F17 (51-100) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	81099 NA 79501 81099 NA<
Gene Symbol	OR4F17 OR4F4 O4F17_HUMAN OR4F4_HUMAN OR4F5_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human OR4F4/4F5/4F17. AA range:51-100
Immunogen Region	51-100
Specificity	Olfactory receptor 4F4/4F5/4F17 Polyclonal Antibody detects endogenous levels of Olfactory receptor 4F4/4F5/4F17 protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081