

Anti-OR4K2 antibody (265-314) (STJ94714)

STJ94714

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Olfactory receptor 4K2 (265-314) is suitable for use in Western Blot, Immunofluorescence and ELISA research applications.
Applications	WB/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	390431
Gene Symbol	OR4K2
Uniprot ID	OR4K2_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human OR4K2. AA range:265-314
Immunogen Region	265-314
Specificity	Olfactory receptor 4K2 Polyclonal Antibody detects endogenous levels of Olfactory receptor 4K2 protein.
Immunogen Sequence	



Western blot analysis of the lysates from HT-29 cells using OR4K2 antibody.



Western blot analysis of lysates from Jurkat cells, using OR4K2 Antibody. The lane on the right is blocked with the synthesized peptide.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081