

Anti-OR52N1 antibody (269-318) (STJ94749)

STJ94749

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Olfactory receptor 52N1 (269-318) is suitable for use in Western Blot, Immunofluorescence and ELISA research applications.
Applications	WB/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION


Gene ID	79473
Gene Symbol	OR52N1
Uniprot ID	O52N1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from human OR52N1. AA range:269-318
Immunogen Region	269-318
Specificity	Olfactory receptor 52N1 Polyclonal Antibody detects endogenous levels of Olfactory receptor 52N1 protein.
Immunogen Sequence	



Western blot analysis of lysates from Jurkat cells, using OR52N1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using OR52N1 antibody.



Immunofluorescence analysis of LOVO cells, using OR52N1 Antibody. The picture on the right is blocked with the synthesized peptide.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081