

Anti-CDKN2D antibody (96-145 aa) (STJ94861)

STJ94861

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cyclin-dependent kinase 4 inhibitor D (96-145 aa) is suitable for use in Immunohistochemistry, Immunofluorescence, Western Blot and ELISA research applications.
Applications	IHC/IF/WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1032
Gene Symbol	CDKN2D
Uniprot ID	CDN2D_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human CDKN2D at the amino acid range 96-145
Immunogen Region	96-145 aa
Specificity	p19 Polyclonal Antibody detects endogenous levels of p19 protein.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human brain, using p19 INK4d Antibody. The picture on the right is blocked with the synthesized peptide.