

Anti-PIAS2 antibody (10-59 aa) (STJ95084)

STJ95084

GENERAL INFORMATION




Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E3 SUMO-protein ligase PIAS2 (10-59 aa) is suitable for use in Western Blot, ELISA and Immunohistochemistry research applications.
Applications	WB/ELISA/IHC
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat/Monkey

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC-P 1:50-300 ELISA 2000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9063
Gene Symbol	PIAS2
Uniprot ID	PIAS2_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human PIAS2 at the amino acid range 10-59
Immunogen Region	10-59 aa
Specificity	PIASx Polyclonal Antibody detects endogenous levels of PIASx protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081