

Anti-PFDN1 antibody (1-50 aa) (STJ95219)

STJ95219

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Prefoldin subunit 1 (1-50 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5201
Gene Symbol	PFDN1
Uniprot ID	PFD1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human PFDN1 at the amino acid range 1-50
Immunogen Region	1-50 aa
Specificity	Prefoldin 1 Polyclonal Antibody detects endogenous levels of Prefoldin 1 protein.
Immunogen Sequence	

