

## Anti-PRDX3 antibody (44-93 aa) (STJ95238)

STJ95238

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Thioredoxin-dependent peroxide reductase, mitochondrial (44-93 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunoprecipitation, Immunofluorescence and ELISA research applications.
<b>Applications</b>	WB/IHC/IP/IF/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES


<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-1:2000
<b>Range</b>	IHC-P 1:50-300 IP 1:50-200 ELISA 1:20000 IF 1:50-200
<b>Formulation</b>	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	10935
<b>Gene Symbol</b>	PRDX3
<b>Uniprot ID</b>	PRDX3_HUMAN
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from the human PRX III at the amino acid range 44-93
<b>Immunogen Region</b>	44-93 aa
<b>Specificity</b>	PRX III Polyclonal Antibody detects endogenous levels of PRX III protein.
<b>Immunogen Sequence</b>	



Immunohistochemistry analysis of PRX III antibody in paraffin-embedded human lung carcinoma tissue.



Western blot analysis of lysate from HeLa cells, using PRX III antibody.



Western blot analysis of various cells using PRX III Polyclonal Antibody diluted at 1/104 1000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081