

Anti-PTGR2 antibody (181-230 aa) (STJ95261)

STJ95261

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Prostaglandin reductase 2 (181-230 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	145482
Gene Symbol	PTGR2
Uniprot ID	PTGR2_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human ZADH1 at the amino acid range 181-230
Immunogen Region	181-230 aa
Specificity	PTGR2 Polyclonal Antibody detects endogenous levels of PTGR2 protein.
Immunogen Sequence	

