

Anti-PTPRA antibody (730-810 aa) (STJ95267)

STJ95267

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Receptor-type tyrosine-protein phosphatase alpha (730-810 aa) is suitable for use in Western Blot, ELISA and Immunohistochemistry research applications.
Applications	WB/ELISA/IHC
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC-P 1:50-300 ELISA 2000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5786
Gene Symbol	PTPRA
Uniprot ID	PTPRA_HUMAN
Immunogen	Synthesized peptide derived from PTP Alpha at the amino acid range 730-810
Immunogen Region	730-810 aa
Specificity	PTP Alpha Polyclonal Antibody detects endogenous levels of PTP Alpha protein.
Immunogen Sequence	



Immunohistochemical analysis of paraffin-embedded human uterus. 1. Antibody was diluted at 1:200 (4A°C overnight). 2. Tris-EDTA, pH9.0 was used for antigen retrieval. 3. Secondary antibody was diluted at 1:200 (room temperature, 45min).