

Anti-Rac1/2/3/CDC42 antibody (38-87 aa) (STJ95319)

STJ95319

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ras-related C3 botulinum toxin substrate 1 and Ras-related C3 botulinum toxin substrate 2 and Ras-related C3 botulinum toxin substrate 3 and Cell division control protein 42 homolog (38-87 aa) is suitable for use in We
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat





PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	IHC 1:100-1:300 ELISA 1:20000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID [5879](#)
[5881](#)
[998](#)
[RAC1](#)
[RAC3](#)
[RAC1_HUMAN](#)
[RAC3_HUMAN](#)
[CDC42_HUMAN](#)

Immunogen The antiserum was produced against synthesized peptide derived from the human Rac1/CDC42 at the amino acid range 38-87
Immunogen Region 38-87 aa
Specificity Rac1/2/3/CDC42 Polyclonal Antibody detects endogenous levels of Rac1/2/3/CDC42 protein.
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081