

Anti-RERG antibody (70-150 aa) (STJ95420)

STJ95420

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ras-related and estrogen-regulated growth inhibitor (70-150 aa) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC 1:100-1:300
Range	ELISA 1:5000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	85004
Gene Symbol	RERG
Uniprot ID	RERG_HUMAN
Immunogen	Synthesized peptide derived from RERG at the amino acid range 70-150
Immunogen Region	70-150 aa
Specificity	RERG Polyclonal Antibody detects endogenous levels of RERG protein.
Immunogen Sequence	